

Package ‘nanoarrow’

December 4, 2025

Title Interface to the 'nanoarrow' 'C' Library

Version 0.7.0-2

Description Provides an 'R' interface to the 'nanoarrow' 'C' library and the 'Apache Arrow' application binary interface. Functions to import and export 'ArrowArray', 'ArrowSchema', and 'ArrowArrayStream' 'C' structures to and from 'R' objects are provided alongside helpers to facilitate zero-copy data transfer among 'R' bindings to libraries implementing the 'Arrow' 'C' data interface.

License Apache License (>= 2)

Encoding UTF-8

RoxygenNote 7.3.2

URL <https://arrow.apache.org/nanoarrow/latest/r/>,
<https://github.com/apache/arrow-nanoarrow>

BugReports <https://github.com/apache/arrow-nanoarrow/issues>

Suggests arrow (>= 9.0.0), bit64, blob, dplyr, hms, jsonlite, rlang,
testthat (>= 3.0.0), tibble, vctrs, withr

SystemRequirements libzstd (optional)

Config/testthat/edition 3

Config/build/bootstrap TRUE

NeedsCompilation yes

Author Dewey Dunnington [aut, cre] (ORCID:
<<https://orcid.org/0000-0002-9415-4582>>),
Apache Arrow [aut, cph],
Apache Software Foundation [cph]

Maintainer Dewey Dunnington <dewey@dunnington.ca>

Repository CRAN

Date/Publication 2025-12-04 15:20:02 UTC

Contents

array_stream_set_finalizer	2
as_nanoarrow_array	3
as_nanoarrow_array_stream	4
as_nanoarrow_buffer	5
as_nanoarrow_schema	5
as_nanoarrow_vctr	6
basic_array_stream	7
convert_array	8
convert_array_stream	10
example_ipc_stream	11
infer_nanoarrow_ptype	11
infer_nanoarrow_ptype_extension	12
nanoarrow_array_init	13
nanoarrow_buffer_init	14
nanoarrow_extension_array	15
nanoarrow_extension_spec	16
nanoarrow_pointer_is_valid	17
nanoarrow_version	19
na_type	19
na_vctrs	23
read_nanoarrow	24
Index	25

array_stream_set_finalizer

Register an array stream finalizer

Description

In some cases, R functions that return a [nanoarrow_array_stream](#) may require that the scope of some other object outlive that of the array stream. If there is a need for that object to be released deterministically (e.g., to close open files), you can register a function to run after the stream's release callback is invoked from the R thread. Note that this finalizer will **not** be run if the stream's release callback is invoked from a **non-R** thread. In this case, the finalizer and its chain of environments will be garbage-collected when `nanoarrow:::preserved_empty()` is run.

Usage

```
array_stream_set_finalizer(array_stream, finalizer)
```

Arguments

`array_stream` A [nanoarrow_array_stream](#)

`finalizer` A function that will be called with zero arguments.

Value

A newly allocated array_stream whose release callback will call the supplied finalizer.

Examples

```
stream <- array_stream_set_finalizer(
  basic_array_stream(list(1:5)),
  function() message("All done!")
)
stream$release()
```

as_nanoarrow_array	<i>Convert an object to a nanoarrow array</i>
--------------------	---

Description

In nanoarrow an 'array' refers to the struct ArrowArray definition in the Arrow C data interface. At the R level, we attach a [schema](#) such that functionally the nanoarrow_array class can be used in a similar way as an arrow::Array. Note that in nanoarrow an arrow::RecordBatch and a non-nullable arrow::StructArray are represented identically.

Usage

```
as_nanoarrow_array(x, ..., schema = NULL)
```

Arguments

x	An object to convert to a array
...	Passed to S3 methods
schema	An optional schema used to enforce conversion to a particular type. Defaults to infer_nanoarrow_schema() .

Value

An object of class 'nanoarrow_array'

Examples

```
(array <- as_nanoarrow_array(1:5))
as.vector(array)

(array <- as_nanoarrow_array(data.frame(x = 1:5)))
as.data.frame(array)
```

`as_nanoarrow_array_stream`*Convert an object to a nanoarrow array_stream*

Description

In nanoarrow, an 'array stream' corresponds to the struct `ArrowArrayStream` as defined in the Arrow C Stream interface. This object is used to represent a stream of [arrays](#) with a common [schema](#). This is similar to an [arrow::RecordBatchReader](#) except it can be used to represent a stream of any type (not just record batches). Note that a stream of record batches and a stream of non-nullable struct arrays are represented identically. Also note that array streams are mutable objects and are passed by reference and not by value.

Usage

```
as_nanoarrow_array_stream(x, ..., schema = NULL)
```

Arguments

<code>x</code>	An object to convert to a <code>array_stream</code>
<code>...</code>	Passed to S3 methods
<code>schema</code>	An optional schema used to enforce conversion to a particular type. Defaults to infer_nanoarrow_schema() .

Value

An object of class 'nanoarrow_array_stream'

Examples

```
(stream <- as_nanoarrow_array_stream(data.frame(x = 1:5)))
stream$get_schema()
stream$get_next()

# The last batch is returned as NULL
stream$get_next()

# Release the stream
stream$release()
```

as_nanoarrow_buffer	<i>Convert an object to a nanoarrow buffer</i>
---------------------	--

Description

Convert an object to a nanoarrow buffer

Usage

```
as_nanoarrow_buffer(x, ...)
```

Arguments

x	An object to convert to a buffer
...	Passed to S3 methods

Value

An object of class 'nanoarrow_buffer'

Examples

```
array <- as_nanoarrow_array(c(NA, 1:4))
array$buffers
as.raw(array$buffers[[1]])
as.raw(array$buffers[[2]])
convert_buffer(array$buffers[[1]])
convert_buffer(array$buffers[[2]])
```

as_nanoarrow_schema	<i>Convert an object to a nanoarrow schema</i>
---------------------	--

Description

In nanoarrow a 'schema' refers to a struct ArrowSchema as defined in the Arrow C Data interface. This data structure can be used to represent an [arrow::schema\(\)](#), an [arrow::field\(\)](#), or an [arrow::DataType](#). Note that in nanoarrow, an [arrow::schema\(\)](#) and a non-nullable [arrow::struct\(\)](#) are represented identically.

Usage

```
as_nanoarrow_schema(x, ...)

infer_nanoarrow_schema(x, ...)

nanoarrow_schema_parse(x, recursive = FALSE)

nanoarrow_schema_modify(x, new_values, validate = TRUE)
```

Arguments

<code>x</code>	An object to convert to a schema
<code>...</code>	Passed to S3 methods
<code>recursive</code>	Use TRUE to include a children member when parsing schemas.
<code>new_values</code>	New schema component to assign
<code>validate</code>	Use FALSE to skip schema validation

Value

An object of class 'nanoarrow_schema'

Examples

```
infer_nanoarrow_schema(integer())
infer_nanoarrow_schema(data.frame(x = integer()))
```

as_nanoarrow_vctr	<i>Experimental Arrow encoded arrays as R vectors</i>
-------------------	---

Description

This experimental vctr class allows zero or more Arrow arrays to present as an R vector without converting them. This is useful for arrays with types that do not have a non-lossy R equivalent, and helps provide an intermediary object type where the default conversion is prohibitively expensive (e.g., a nested list of data frames). These objects will not survive many vctr transformations; however, they can be sliced without copying the underlying arrays.

Usage

```
as_nanoarrow_vctr(x, ..., schema = NULL, subclass = character())

nanoarrow_vctr(schema = NULL, subclass = character())
```

Arguments

x	An object that works with <code>as_nanoarrow_array_stream()</code> .
...	Passed to <code>as_nanoarrow_array_stream()</code>
schema	An optional schema
subclass	An optional subclass of <code>nanoarrow_vctr</code> to prepend to the final class name.

Details

The `nanoarrow_vctr` is currently implemented similarly to `factor()`: its storage type is an `integer()` that is a sequence along the total length of the `vctr` and there are attributes that are required to resolve these indices to an array + offset. Sequences typically have a very compact representation in recent versions of R such that this has a cheap storage footprint even for large arrays. The attributes are currently:

- schema: The `nanoarrow_schema` shared by each chunk.
- chunks: A `list()` of `nanoarrow_array`.
- offsets: An `integer()` vector beginning with 0 and followed by the cumulative length of each chunk. This allows the chunk index + offset to be resolved from a logical index with $\log(n)$ complexity.

This implementation is preliminary and may change; however, the result of `as_nanoarrow_array_stream(some_vctr[begin:end])` should remain stable.

Value

A `vctr` of class `'nanoarrow_vctr'`

Examples

```
array <- as_nanoarrow_array(1:5)
as_nanoarrow_vctr(array)
```

<code>basic_array_stream</code>	<i>Create ArrayStreams from batches</i>
---------------------------------	---

Description

Create ArrayStreams from batches

Usage

```
basic_array_stream(batches, schema = NULL, validate = TRUE)
```

Arguments

batches	A <code>list()</code> of <code>nanoarrow_array</code> objects or objects that can be coerced via <code>as_nanoarrow_array()</code> .
schema	A <code>nanoarrow_schema</code> or <code>NULL</code> to guess based on the first schema.
validate	Use <code>FALSE</code> to skip the validation step (i.e., if you know that the arrays are valid).

Value

An `nanoarrow_array_stream`

Examples

```
(stream <- basic_array_stream(list(data.frame(a = 1, b = 2)))
as.data.frame(stream$get_next())
stream$get_next())
```

convert_array

Convert an Array into an R vector

Description

Converts array to the type specified by `to`. This is a low-level interface; most users should use `as.data.frame()` or `as.vector()` unless finer-grained control is needed over the conversion. This function is an S3 generic dispatching on `to`: developers may implement their own S3 methods for custom vector types.

Usage

```
convert_array(array, to = NULL, ...)
```

Arguments

array	A <code>nanoarrow_array</code> .
to	A target prototype object describing the type to which array should be converted, or <code>NULL</code> to use the default conversion as returned by <code>infer_nanoarrow_ptype()</code> . Alternatively, a function can be passed to perform an alternative calculation of the default ptype as a function of array and the default inference of the prototype.
...	Passed to S3 methods

Details

Note that unregistered extension types will by default issue a warning. Use `options(nanoarrow.warn_unregistered_extensions = FALSE)` to disable this behaviour.

Conversions are implemented for the following R vector types:

- `logical()`: Any numeric type can be converted to `logical()` in addition to the bool type. For numeric types, any non-zero value is considered TRUE.
- `integer()`: Any numeric type can be converted to `integer()`; however, a warning will be signaled if the any value is outside the range of the 32-bit integer.
- `double()`: Any numeric type can be converted to `double()`. This conversion currently does not warn for values that may not roundtrip through a floating-point double (e.g., very large uint64 and int64 values).
- `character()`: String and large string types can be converted to `character()`. The conversion does not check for valid UTF-8: if you need finer-grained control over encodings, use `to = blob::blob()`.
- `factor()`: Dictionary-encoded arrays of strings can be converted to `factor()`; however, this must be specified explicitly (i.e., `convert_array(array, factor())`) because arrays arriving in chunks can have dictionaries that contain different levels. Use `convert_array(array, factor(levels = c(...)))` to materialize an array into a vector with known levels.
- `Date`: Only the date32 type can be converted to an R Date vector.
- `hms::hms()`: Time32 and time64 types can be converted to `hms::hms()`.
- `difftime()`: Time32, time64, and duration types can be converted to R `difftime()` vectors. The value is converted to match the `units()` attribute of `to`.
- `blob::blob()`: String, large string, binary, and large binary types can be converted to `blob::blob()`.
- `vctrs::list_of()`: List, large list, and fixed-size list types can be converted to `vctrs::list_of()`.
- `matrix()`: Fixed-size list types can be converted to `matrix(ptype, ncol = fixed_size)`.
- `data.frame()`: Struct types can be converted to `data.frame()`.
- `vctrs::unspecified()`: Any type can be converted to `vctrs::unspecified()`; however, a warning will be raised if any non-null values are encountered.

In addition to the above conversions, a null array may be converted to any target prototype except `data.frame()`. Extension arrays are currently converted as their storage type.

Value

An R vector of type `to`.

Examples

```
array <- as_nanoarrow_array(data.frame(x = 1:5))
str(convert_array(array))
str(convert_array(array, to = data.frame(x = double())))
```

convert_array_stream *Convert an Array Stream into an R vector*

Description

Converts array_stream to the type specified by to. This is a low-level interface; most users should use `as.data.frame()` or `as.vector()` unless finer-grained control is needed over the conversion. See `convert_array()` for details of the conversion process; see `infer_nanoarrow_ptype()` for default inferences of to.

Usage

```
convert_array_stream(array_stream, to = NULL, size = NULL, n = Inf)
```

```
collect_array_stream(array_stream, n = Inf, schema = NULL, validate = TRUE)
```

Arguments

array_stream	A nanoarrow_array_stream .
to	A target prototype object describing the type to which array should be converted, or NULL to use the default conversion as returned by <code>infer_nanoarrow_ptype()</code> . Alternatively, a function can be passed to perform an alternative calculation of the default ptype as a function of array and the default inference of the prototype.
size	The exact size of the output, if known. If specified, slightly more efficient implementation may be used to collect the output.
n	The maximum number of batches to pull from the array stream.
schema	A nanoarrow_schema or NULL to guess based on the first schema.
validate	Use FALSE to skip the validation step (i.e., if you know that the arrays are valid).

Value

- `convert_array_stream()`: An R vector of type to.
- `collect_array_stream()`: A `list()` of [nanoarrow_array](#)

Examples

```
stream <- as_nanoarrow_array_stream(data.frame(x = 1:5))
str(convert_array_stream(stream))
str(convert_array_stream(stream, to = data.frame(x = double()))))

stream <- as_nanoarrow_array_stream(data.frame(x = 1:5))
collect_array_stream(stream)
```

example_ipc_stream	<i>Example Arrow IPC Data</i>
--------------------	-------------------------------

Description

An example stream that can be used for testing or examples.

Usage

```
example_ipc_stream(compression = c("none", "zstd"))
```

Arguments

compression One of "none" or "zstd"

Value

A raw vector that can be passed to [read_nanoarrow\(\)](#)

Examples

```
as.data.frame(read_nanoarrow(example_ipc_stream()))
```

infer_nanoarrow_ptype	<i>Infer an R vector prototype</i>
-----------------------	------------------------------------

Description

Resolves the default to value to use in [convert_array\(\)](#) and [convert_array_stream\(\)](#). The default conversions are:

Usage

```
infer_nanoarrow_ptype(x)
```

Arguments

x A [nanoarrow_schema](#), [nanoarrow_array](#), or [nanoarrow_array_stream](#).

Details

- null to `vctrs::unspecified()`
- boolean to `logical()`
- int8, uint8, int16, uint16, and int32 to `integer()`
- uint32, int64, uint64, float, and double to `double()`
- string and large string to `character()`
- struct to `data.frame()`
- binary and large binary to `blob::blob()`
- list, large_list, and fixed_size_list to `vctrs::list_of()`
- time32 and time64 to `hms::hms()`
- duration to `difftime()`
- date32 to `as.Date()`
- timestamp to `as.POSIXct()`

Additional conversions are possible by specifying an explicit value for `to`. For details of each conversion, see `convert_array()`.

Value

An R vector of zero size describing the target into which the array should be materialized.

Examples

```
infer_nanoarrow_ptype(as_nanoarrow_array(1:10))
```

```
infer_nanoarrow_ptype_extension
```

Implement Arrow extension types

Description

Implement Arrow extension types

Usage

```
infer_nanoarrow_ptype_extension(
  extension_spec,
  x,
  ...,
  warn_unregistered = TRUE
)

convert_array_extension(
```

```

    extension_spec,
    array,
    to,
    ...,
    warn_unregistered = TRUE
)

as_nanoarrow_array_extension(extension_spec, x, ..., schema = NULL)

```

Arguments

`extension_spec` An extension specification inheriting from 'nanoarrow_extension_spec'.

`x, array, to, schema, ...`
 Passed from [infer_nanoarrow_ptype\(\)](#), [convert_array\(\)](#), [as_nanoarrow_array\(\)](#),
 and/or [as_nanoarrow_array_stream\(\)](#).

`warn_unregistered`
 Use FALSE to infer/convert based on the storage type without a warning.

Value

- `infer_nanoarrow_ptype_extension()`: The R vector prototype to be used as the default conversion target.
- `convert_array_extension()`: An R vector of type `to`.
- `as_nanoarrow_array_extension()`: A [nanoarrow_array](#) of type `schema`.

`nanoarrow_array_init` *Modify nanoarrow arrays*

Description

Create a new array or from an existing array, modify one or more parameters. When importing an array from elsewhere, `nanoarrow_array_set_schema()` is useful to attach the data type information to the array (without this information there is little that nanoarrow can do with the array since its content cannot be otherwise interpreted). `nanoarrow_array_modify()` can create a shallow copy and modify various parameters to create a new array, including setting children and buffers recursively. These functions power the `$<-` operator, which can modify one parameter at a time.

Usage

```

nanoarrow_array_init(schema)

nanoarrow_array_set_schema(array, schema, validate = TRUE)

nanoarrow_array_modify(array, new_values, validate = TRUE)

```

Arguments

schema	A nanoarrow_schema to attach to this array.
array	A nanoarrow_array .
validate	Use FALSE to skip validation. Skipping validation may result in creating an array that will crash R.
new_values	A named <code>list()</code> of values to replace.

Value

- `nanoarrow_array_init()` returns a possibly invalid but initialized array with a given schema.
- `nanoarrow_array_set_schema()` returns array, invisibly. Note that array is modified in place by reference.
- `nanoarrow_array_modify()` returns a shallow copy of array with the modified parameters such that the original array remains valid.

Examples

```
nanoarrow_array_init(na_string())

# Modify an array using $ and <-
array <- as_nanoarrow_array(1:5)
array$length <- 4
as.vector(array)

# Modify potentially more than one component at a time
array <- as_nanoarrow_array(1:5)
as.vector(nanoarrow_array_modify(array, list(length = 4)))

# Attach a schema to an array
array <- as_nanoarrow_array(-1L)
nanoarrow_array_set_schema(array, na_uint32())
as.vector(array)
```

`nanoarrow_buffer_init` *Create and modify nanoarrow buffers*

Description

Create and modify nanoarrow buffers

Usage

```
nanoarrow_buffer_init()

nanoarrow_buffer_append(buffer, new_buffer)

convert_buffer(buffer, to = NULL)
```

Arguments

buffer, new_buffer
[nanoarrow_buffers](#).

to A target prototype object describing the type to which array should be converted, or NULL to use the default conversion as returned by [infer_nanoarrow_ptype\(\)](#). Alternatively, a function can be passed to perform an alternative calculation of the default ptype as a function of array and the default inference of the prototype.

Value

- `nanoarrow_buffer_init()`: An object of class 'nanoarrow_buffer'
- `nanoarrow_buffer_append()`: Returns buffer, invisibly. Note that buffer is modified in place by reference.

Examples

```
buffer <- nanoarrow_buffer_init()
nanoarrow_buffer_append(buffer, 1:5)

array <- nanoarrow_array_modify(
  nanoarrow_array_init(na_int32()),
  list(length = 5, buffers = list(NULL, buffer))
)
as.vector(array)
```

nanoarrow_extension_array

Create Arrow extension arrays

Description

Create Arrow extension arrays

Usage

```
nanoarrow_extension_array(
  storage_array,
  extension_name,
  extension_metadata = NULL
)
```

Arguments

storage_array A [nanoarrow_array](#).

extension_name For [na_extension\(\)](#), the extension name. This is typically namespaced separated by dots (e.g., nanoarrow.r.vctrs).

extension_metadata
A string or raw vector defining extension metadata. Most Arrow extension types define extension metadata as a JSON object.

Value

A [nanoarrow_array](#) with attached extension schema.

Examples

```
nanoarrow_extension_array(1:10, "some_ext", '{"key": "value"}')
```

nanoarrow_extension_spec

Register Arrow extension types

Description

Register Arrow extension types

Usage

```
nanoarrow_extension_spec(data = list(), subclass = character())

register_nanoarrow_extension(extension_name, extension_spec)

unregister_nanoarrow_extension(extension_name)

resolve_nanoarrow_extension(extension_name)
```

Arguments

data Optional data to include in the extension type specification

subclass A subclass for the extension type specification. Extension methods will dispatch on this object.

extension_name An Arrow extension type name (e.g., nanoarrow.r.vctrs)

extension_spec An extension specification inheriting from 'nanoarrow_extension_spec'.

Value

- `nanarrow_extension_spec()` returns an object of class 'nanarrow_extension_spec'.
- `register_nanarrow_extension()` returns `extension_spec`, invisibly.
- `unregister_nanarrow_extension()` returns `extension_name`, invisibly.
- `resolve_nanarrow_extension()` returns an object of class 'nanarrow_extension_spec' or NULL if the extension type was not registered.

Examples

```
nanarrow_extension_spec("mynamespace.mytype", subclass = "mypackage_mytype_spec")
```

```
nanarrow_pointer_is_valid
```

Danger zone: low-level pointer operations

Description

The [nanarrow_schema](#), [nanarrow_array](#), and [nanarrow_array_stream](#) classes are represented in R as external pointers (EXTPTRXP). When these objects go out of scope (i.e., when they are garbage collected or shortly thereafter), the underlying object's `release()` callback is called if the underlying pointer is non-null and if the `release()` callback is non-null.

Usage

```
nanarrow_pointer_is_valid(ptr)

nanarrow_pointer_addr_dbl(ptr)

nanarrow_pointer_addr_chr(ptr)

nanarrow_pointer_addr_pretty(ptr)

nanarrow_pointer_release(ptr)

nanarrow_pointer_move(ptr_src, ptr_dst)

nanarrow_pointer_export(ptr_src, ptr_dst)

nanarrow_allocate_schema()

nanarrow_allocate_array()

nanarrow_allocate_array_stream()

nanarrow_pointer_set_protected(ptr_src, protected)
```

Arguments

<code>ptr, ptr_src, ptr_dst</code>	An external pointer to a struct ArrowSchema, struct ArrowArray, or struct ArrowArrayStream.
<code>protected</code>	An object whose scope must outlive that of <code>ptr</code> . This is useful for array streams since at least two specifications involving the array stream specify that the stream is only valid for the lifecycle of another object (e.g., an AdbcStatement or OGR-Dataset).

Details

When interacting with other C Data Interface implementations, it is important to keep in mind that the R object wrapping these pointers is always passed by reference (because it is an external pointer) and may be referred to by another R object (e.g., an element in a `list()` or as a variable assigned in a user's environment). When importing a schema, array, or array stream into nanoarrow this is not a problem: the R object takes ownership of the lifecycle and memory is released when the R object is garbage collected. In this case, one can use `nanoarrow_pointer_move()` where `ptr_dst` was created using `nanoarrow_allocate_*`.

The case of exporting is more complicated and as such has a dedicated function, `nanoarrow_pointer_export()`, that implements different logic schemas, arrays, and array streams:

- Schema objects are (deep) copied such that a fresh copy of the schema is exported and made the responsibility of some other C data interface implementation.
- Array objects are exported as a shell around the original array that preserves a reference to the R object. This ensures that the buffers and children pointed to by the array are not copied and that any references to the original array are not invalidated.
- Array stream objects are moved: the responsibility for the object is transferred to the other C data interface implementation and any references to the original R object are invalidated. Because these objects are mutable, this is typically what you want (i.e., you should not be pulling arrays from a stream accidentally from two places).

If you know the lifecycle of your object (i.e., you created the R object yourself and never passed references to it elsewhere), you can slightly more efficiently call `nanoarrow_pointer_move()` for all three pointer types.

Value

- `nanoarrow_pointer_is_valid()` returns TRUE if the pointer is non-null and has a non-null release callback.
- `nanoarrow_pointer_addr_dbl()` and `nanoarrow_pointer_addr_chr()` return pointer representations that may be helpful to facilitate moving or exporting nanoarrow objects to other libraries.
- `nanoarrow_pointer_addr_pretty()` gives a pointer representation suitable for printing or error messages.
- `nanoarrow_pointer_release()` returns `ptr`, invisibly.
- `nanoarrow_pointer_move()` and `nanoarrow_pointer_export()` return `ptr_dst`, invisibly.

- nanoarrow_allocate_array(), nanoarrow_allocate_schema(), and nanoarrow_allocate_array_stream() return an [array](#), a [schema](#), and an [array stream](#), respectively.

nanoarrow_version	<i>Underlying 'nanoarrow' C library build</i>
-------------------	---

Description

Underlying 'nanoarrow' C library build

Usage

```
nanoarrow_version(runtime = TRUE)

nanoarrow_with_zstd()
```

Arguments

runtime Compare TRUE and FALSE values to detect a possible ABI mismatch.

Value

A string identifying the version of nanoarrow this package was compiled against.

Examples

```
nanoarrow_version()
nanoarrow_with_zstd()
```

na_type	<i>Create type objects</i>
---------	----------------------------

Description

In nanoarrow, types, fields, and schemas are all represented by a [nanoarrow_schema](#). These functions are convenience constructors to create these objects in a readable way. Use [na_type\(\)](#) to construct types based on the constructor name, which is also the name that prints/is returned by [nanoarrow_schema_parse\(\)](#).

Usage

```
na_type(  
  type_name,  
  byte_width = NULL,  
  unit = NULL,  
  timezone = NULL,  
  precision = NULL,  
  scale = NULL,  
  column_types = NULL,  
  item_type = NULL,  
  key_type = NULL,  
  value_type = NULL,  
  index_type = NULL,  
  ordered = NULL,  
  list_size = NULL,  
  keys_sorted = NULL,  
  storage_type = NULL,  
  extension_name = NULL,  
  extension_metadata = NULL,  
  nullable = NULL  
)  
  
na_na(nullable = TRUE)  
  
na_bool(nullable = TRUE)  
  
na_int8(nullable = TRUE)  
  
na_uint8(nullable = TRUE)  
  
na_int16(nullable = TRUE)  
  
na_uint16(nullable = TRUE)  
  
na_int32(nullable = TRUE)  
  
na_uint32(nullable = TRUE)  
  
na_int64(nullable = TRUE)  
  
na_uint64(nullable = TRUE)  
  
na_half_float(nullable = TRUE)  
  
na_float(nullable = TRUE)  
  
na_double(nullable = TRUE)
```

```
na_string(nullable = TRUE)
na_large_string(nullable = TRUE)
na_string_view(nullable = TRUE)
na_binary(nullable = TRUE)
na_large_binary(nullable = TRUE)
na_fixed_size_binary(byte_width, nullable = TRUE)
na_binary_view(nullable = TRUE)
na_date32(nullable = TRUE)
na_date64(nullable = TRUE)
na_time32(unit = c("ms", "s"), nullable = TRUE)
na_time64(unit = c("us", "ns"), nullable = TRUE)
na_duration(unit = c("ms", "s", "us", "ns"), nullable = TRUE)
na_interval_months(nullable = TRUE)
na_interval_day_time(nullable = TRUE)
na_interval_month_day_nano(nullable = TRUE)
na_timestamp(unit = c("us", "ns", "s", "ms"), timezone = "", nullable = TRUE)
na_decimal32(precision, scale, nullable = TRUE)
na_decimal64(precision, scale, nullable = TRUE)
na_decimal128(precision, scale, nullable = TRUE)
na_decimal256(precision, scale, nullable = TRUE)
na_struct(column_types = list(), nullable = FALSE)
na_sparse_union(column_types = list())
na_dense_union(column_types = list())
na_list(item_type, nullable = TRUE)
```

```

na_large_list(item_type, nullable = TRUE)

na_list_view(item_type, nullable = TRUE)

na_large_list_view(item_type, nullable = TRUE)

na_fixed_size_list(item_type, list_size, nullable = TRUE)

na_map(key_type, item_type, keys_sorted = FALSE, nullable = TRUE)

na_dictionary(value_type, index_type = na_int32(), ordered = FALSE)

na_extension(storage_type, extension_name, extension_metadata = "")

```

Arguments

type_name	The name of the type (e.g., "int32"). This form of the constructor is useful for writing tests that loop over many types.
byte_width	For na_fixed_size_binary() , the number of bytes occupied by each item.
unit	One of 's' (seconds), 'ms' (milliseconds), 'us' (microseconds), or 'ns' (nanoseconds).
timezone	A string representing a timezone name. The empty string "" represents a naive point in time (i.e., one that has no associated timezone).
precision	The total number of digits representable by the decimal type
scale	The number of digits after the decimal point in a decimal type
column_types	A list() of nanoarrow_schemas .
item_type	For na_list() , na_large_list() , na_fixed_size_list() , and na_map() , the nanoarrow_schema representing the item type.
key_type	The nanoarrow_schema representing the na_map() key type.
value_type	The nanoarrow_schema representing the na_dictionary() or na_map() value type.
index_type	The nanoarrow_schema representing the na_dictionary() index type.
ordered	Use TRUE to assert that the order of values in the dictionary are meaningful.
list_size	The number of elements in each item in a na_fixed_size_list() .
keys_sorted	Use TRUE to assert that keys are sorted.
storage_type	For na_extension() , the underlying value type.
extension_name	For na_extension() , the extension name. This is typically namespaced separated by dots (e.g., nanoarrow.r.vctrs).
extension_metadata	A string or raw vector defining extension metadata. Most Arrow extension types define extension metadata as a JSON object.
nullable	Use FALSE to assert that this field cannot contain null values.

Value

A [nanoarrow_schema](#)

Examples

```
na_int32()
na_struct(list(col1 = na_int32()))
```

na_vctrs	<i>Vctrs extension type</i>
----------	-----------------------------

Description

The Arrow format provides a rich type system that can handle most R vector types; however, many R vector types do not roundtrip perfectly through Arrow memory. The vctrs extension type uses [vctrs::vec_data\(\)](#), [vctrs::vec_restore\(\)](#), and [vctrs::vec_ptype\(\)](#) in calls to [as_nanoarrow_array\(\)](#) and [convert_array\(\)](#) to ensure roundtrip fidelity.

Usage

```
na_vctrs(ptype, storage_type = NULL)
```

Arguments

ptype	A vctrs prototype as returned by vctrs::vec_ptype() . The prototype can be of arbitrary size, but a zero-size vector is sufficient here.
storage_type	For na_extension() , the underlying value type.

Value

A [nanoarrow_schema](#).

Examples

```
vctr <- as.POSIXlt("2000-01-02 03:45", tz = "UTC")
array <- as_nanoarrow_array(vctr, schema = na_vctrs(vctr))
infer_nanoarrow_ptype(array)
convert_array(array)
```

read_nanoarrow

*Read/write serialized streams of Arrow data***Description**

Reads/writes connections, file paths, URLs, or raw vectors from/to serialized Arrow data. Arrow documentation typically refers to this format as "Arrow IPC", since its origin was as a means to transmit tables between processes (e.g., multiple R sessions). This format can also be written to and read from files or URLs and is essentially a high performance equivalent of a CSV file that does a better job maintaining types.

Usage

```
read_nanoarrow(x, ..., lazy = FALSE)
```

```
write_nanoarrow(data, x, ...)
```

Arguments

<code>x</code>	A <code>raw()</code> vector, connection, or file path from which to read binary data. Common extensions indicating compression (<code>.gz</code> , <code>.bz2</code> , <code>.zip</code>) are automatically uncompressed.
<code>...</code>	Currently unused.
<code>lazy</code>	By default, <code>read_nanoarrow()</code> will read and discard a copy of the reader's schema to ensure that invalid streams are discovered as soon as possible. Use <code>lazy = TRUE</code> to defer this check until the reader is actually consumed.
<code>data</code>	An object to write as an Arrow IPC stream, converted using <code>as_nanoarrow_array_stream()</code> . Notably, this includes a <code>data.frame()</code> .

Details

The nanoarrow package implements an IPC writer; however, you can also use `arrow::write_ipc_stream()` to write data from R, or use the equivalent writer from another Arrow implementation in Python, C++, Rust, JavaScript, Julia, C#, and beyond.

The media type of an Arrow stream is `application/vnd.apache.arrow.stream` and the recommended file extension is `.arrows`.

Value

A `nanoarrow_array_stream`

Examples

```
as.data.frame(read_nanoarrow(example_ipc_stream()))
```


Index

array, [19](#)
array_stream, [19](#)
array_stream_set_finalizer, [2](#)
arrays, [4](#)
arrow::field(), [5](#)
arrow::RecordBatchReader, [4](#)
arrow::schema(), [5](#)
arrow::struct(), [5](#)
arrow::write_ipc_stream(), [24](#)
as.Date(), [12](#)
as.POSIXct(), [12](#)
as_nanoarrow_array, [3](#)
as_nanoarrow_array(), [8](#), [13](#), [23](#)
as_nanoarrow_array_extension
 (infer_nanoarrow_ptype_extension),
 [12](#)
as_nanoarrow_array_stream, [4](#)
as_nanoarrow_array_stream(), [7](#), [13](#), [24](#)
as_nanoarrow_buffer, [5](#)
as_nanoarrow_schema, [5](#)
as_nanoarrow_vctr, [6](#)

basic_array_stream, [7](#)
blob::blob(), [9](#), [12](#)

character(), [9](#), [12](#)
collect_array_stream
 (convert_array_stream), [10](#)
convert_array, [8](#)
convert_array(), [10–13](#), [23](#)
convert_array_extension
 (infer_nanoarrow_ptype_extension),
 [12](#)
convert_array_stream, [10](#)
convert_array_stream(), [11](#)
convert_buffer(nanoarrow_buffer_init),
 [14](#)

data.frame(), [9](#), [12](#), [24](#)
Date, [9](#)

difftime(), [9](#), [12](#)
double(), [9](#), [12](#)

example_ipc_stream, [11](#)

factor(), [9](#)

hms::hms(), [9](#), [12](#)

infer_nanoarrow_ptype, [11](#)
infer_nanoarrow_ptype(), [8](#), [10](#), [13](#), [15](#)
infer_nanoarrow_ptype_extension, [12](#)
infer_nanoarrow_schema
 (as_nanoarrow_schema), [5](#)
infer_nanoarrow_schema(), [3](#), [4](#)
integer(), [9](#), [12](#)

list(), [8](#)
logical(), [9](#), [12](#)

matrix(), [9](#)

na_binary(na_type), [19](#)
na_binary_view(na_type), [19](#)
na_bool(na_type), [19](#)
na_date32(na_type), [19](#)
na_date64(na_type), [19](#)
na_decimal128(na_type), [19](#)
na_decimal256(na_type), [19](#)
na_decimal32(na_type), [19](#)
na_decimal64(na_type), [19](#)
na_dense_union(na_type), [19](#)
na_dictionary(na_type), [19](#)
na_dictionary(), [22](#)
na_double(na_type), [19](#)
na_duration(na_type), [19](#)
na_extension(na_type), [19](#)
na_extension(), [16](#), [22](#), [23](#)
na_fixed_size_binary(na_type), [19](#)
na_fixed_size_binary(), [22](#)
na_fixed_size_list(na_type), [19](#)

- na_fixed_size_list(), 22
- na_float (na_type), 19
- na_half_float (na_type), 19
- na_int16 (na_type), 19
- na_int32 (na_type), 19
- na_int64 (na_type), 19
- na_int8 (na_type), 19
- na_interval_day_time (na_type), 19
- na_interval_month_day_nano (na_type), 19
- na_interval_months (na_type), 19
- na_large_binary (na_type), 19
- na_large_list (na_type), 19
- na_large_list(), 22
- na_large_list_view (na_type), 19
- na_large_string (na_type), 19
- na_list (na_type), 19
- na_list(), 22
- na_list_view (na_type), 19
- na_map (na_type), 19
- na_map(), 22
- na_na (na_type), 19
- na_sparse_union (na_type), 19
- na_string (na_type), 19
- na_string_view (na_type), 19
- na_struct (na_type), 19
- na_time32 (na_type), 19
- na_time64 (na_type), 19
- na_timestamp (na_type), 19
- na_type, 19
- na_type(), 19
- na_uint16 (na_type), 19
- na_uint32 (na_type), 19
- na_uint64 (na_type), 19
- na_uint8 (na_type), 19
- na_vctrs, 23
- nanoarrow_allocate_array
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_allocate_array_stream
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_allocate_schema
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_array, 8, 10, 11, 13, 14, 16, 17
- nanoarrow_array_init, 13
- nanoarrow_array_modify
 - (nanoarrow_array_init), 13
- nanoarrow_array_set_schema
 - (nanoarrow_array_init), 13
- nanoarrow_array_stream, 2, 8, 10, 11, 17, 24
- nanoarrow_buffer, 15
- nanoarrow_buffer_append
 - (nanoarrow_buffer_init), 14
- nanoarrow_buffer_init, 14
- nanoarrow_extension_array, 15
- nanoarrow_extension_spec, 16
- nanoarrow_pointer_addr_chr
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_addr_dbl
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_addr_pretty
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_export
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_export(), 18
- nanoarrow_pointer_is_valid, 17
- nanoarrow_pointer_move
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_move(), 18
- nanoarrow_pointer_release
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_pointer_set_protected
 - (nanoarrow_pointer_is_valid), 17
- nanoarrow_schema, 7, 8, 10, 11, 14, 17, 19, 22, 23
- nanoarrow_schema_modify
 - (as_nanoarrow_schema), 5
- nanoarrow_schema_parse
 - (as_nanoarrow_schema), 5
- nanoarrow_schema_parse(), 19
- nanoarrow_vctr (as_nanoarrow_vctr), 6
- nanoarrow_version, 19
- nanoarrow_with_zstd
 - (nanoarrow_version), 19
- read_nanoarrow, 24
- read_nanoarrow(), 11
- register_nanoarrow_extension
 - (nanoarrow_extension_spec), 16

`resolve_nanoarrow_extension`
 (`nanoarrow_extension_spec`), 16

`schema`, 3, 4, 19

`units()`, 9

`unregister_nanoarrow_extension`
 (`nanoarrow_extension_spec`), 16

`vctrs::list_of()`, 9, 12

`vctrs::unspecified()`, 9, 12

`vctrs::vec_data()`, 23

`vctrs::vec_ptype()`, 23

`vctrs::vec_restore()`, 23

`write_nanoarrow(read_nanoarrow)`, 24